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Moorfields
Eye Hospitals UAE
A branch of Moorfields London

Schwind Amaris 750S

When you decide on an eye laser treatment, you expect the best possible results. The more fully developed the methods are, the better the outcome will be. The SCHWIND AMARIS 750S offers you the leading technology for your laser treatment – superior in all important aspects: Speed, precision, safety and comfort.

Minimum treatment time

The SCHWIND AMARIS 750S works two to three times faster than other laser systems. 750 tiny light flashes per second shape the corneal surface. It is extremely fast and corrects one dioptre of myopia within 1.5 seconds, and eight dioptres are ablated within 13 seconds.

Automatic Fluence Level adjustment

In order to perfectly combine speed and precision, the SCHWIND AMARIS 750S is equipped – besides a 750 Hz pulse frequency – with two energy levels.

At first, the SCHWIND AMARIS 750S uses a high energy level and rapidly ablates around 80% of the tissue to be removed – like a fast yacht quickly heading for the safe harbour.

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

How long the powerful beams are applied depends on your individual refractive error and is continuously checked by the software during the treatment. Like a small sport boat then exactly navigates to a narrow landing, a gentler beam ablates the remaining 20%, thus achieving the especially smooth surface required for excellent vision.

Eye tracking in six dimensions

The eyes can involuntarily move in all directions for milliseconds while fixating on the laser light. The SCHWIND AMARIS 750S eye tracker monitors the position of the eye with approximately 1050 measurements per second and thereby compensates any deviation up to the sixth dimension.

Imagine your eyeball being a buoy in the sea: when the sea is calm, the buoy shifts from the left to the right or forward and backward. In swell, however, the buoy tilts to the side. Furthermore, a buoy also turns around its anchorage or moves upward and downward.

In the same manner, your eyeball can roll (rotation), rotate around its own axis (cyclotorsion) or move up and down along the z-axis. The eye tracker of the SCHWIND AMARIS 750S detects all these eye movements and compensates for them. With cyclotorsion control, the compensation occurs statically between the sitting and the lying positions of the patient, and dynamically during the laser treatment itself.

Permanent measurement of corneal thickness


Similar to a diver orientating himself in water with his depth gauge, your doctor is permanently informed through the online pachymetry of the SCHWIND AMARIS 750S.

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
This system measures the changes of corneal thickness in real time and displays them on the treatment screen. Thus your doctor knows at all times how much tissue has been ablated and the remaining thickness of the cornea.

Exceptionally fine laser beam

The laser beams of the SCHWIND AMARIS 750S have an exceptionally small diameter of only 0.54 millimeters and are thus significantly smaller than others. Furthermore they have a special shape. Your cornea becomes especially smooth with this fine beam.

Constant microclimate

With the SCHWIND AMARIS 750S, two specially formed nozzles ensure a directed, fine airflow that creates a stable microclimate across the cornea. The climatic conditions remain constant, preventing the cornea from drying out.

Unique thermal control

The SCHWIND AMARIS 750S is equipped with a unique thermal control system. Even at high ablation speed the corneal tissue is effectively preserved. The laser pulses are temporally and spatially sorted, which gives the individual position on the cornea sufficient time to cool off. Following pulses can thereby approach an already cooled position faster.

Clinical results

A multitude of international clinical studies have documented the excellent treatment quality with SCHWIND AMARIS technology. Visual acuity of 6/6 (100%) or better was achieved in the vast majority of the cases. A high percentage of treated patients could see better than what they saw before with glasses or with contact lenses. The study results also showed improved contrast vision.

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www.schwindamaris.com

Technology at its best

750 Hz pulse frequency / Automatic Fluence Level Adjustment

▶ Maximum precision, minimum treatment time, 1050 Hz and 6D eye tracker

▶ Registration and compensation of the slightest eye movements up to the sixth dimension

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▶ Gentle ablation of corneal tissue

www.schwindamaris.com

**Moorfields
Eye Hospital Dubai**
A branch of Moorfields London

Moorfields Eye Hospital Dubai

Dubai Healthcare City,
PO Box 505054, District 1,
Al Razi Building 64, Block E, Floor 3,
Dubai, U.A.E.
Tel. +971 4 429 7888
moorfields.ae

Opening hours:

Saturday to Thursday, 8.30am to 5.30pm,
for information and advice on eye conditions and
treatments from experienced ophthalmic trained staff.

Moorfields Abu Dhabi

**Moorfields
Eye Hospital Centre**
A branch of Moorfields London

IN PARTNERSHIP WITH UEMEDICAL

Moorfields Eye Hospital Centre

PO Box 62807, Villas B01 and B02,
Marina Village, near Marina Mall,
Abu Dhabi, U.A.E.
Tel. +971 2 635 6161
moorfields.ae

Opening hours:

Sunday to Wednesday, 9.00am to 8.00pm.
Thursday and Saturday, 9.00am to 6.00pm.

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